


## Spare parts 18510000 Version A



Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza  
Diseños para servicios de reparaciones • Загрузка чертежей для ремонтных мастерских

1851 X000 Easyclean MD



21-1286 24102024/A

Position	Product name	Item number
1	Locking mechanism 1851	900033827
2	Mains cable 230 V with low-power plug	900032750
101	Circuit board configured Easyclean MD	900033832
102	Holder for printed circuit board	900033833
103	Electronics circuit board complete	900033830
104	Ultrasound well 230 V	900033826
105	Seal 1851	900033829
106	Cable harness Easyclean MD	900033831
107	Fuse 5x20mm, T15A	900033837

Position	Product name	Item number
108	Equipment base	900033836
501	Plastic lid	18500001
502	Stainless steel insert cover	18500002
503	Stainless steel basket	18500003
504	Stainless steel immersion basket 59mm	18500004
505	Plastic bath insert for acids	18500005
506	Cleaning jar	18500006
507	Plastic container with lid	18500007 *

\* Wear and tear parts, excluded from the guarantee.