


## Spare parts 18520000 Version A



Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza  
Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

1852 \*000 Easyclean TEC



20-8128 2311/2023 / A

Position	Product name	Item number
1	Rotary knob	900033889
2	Front foil 1852	900033894
3A	Mains cable 230 V with low-power plug	900032750
3B	Mains cable 120V with low-power plug	900035570
101A	Circuit card Easyclean TEC 220-240 V	900033887
101B	Circuit card Easyclean TEC 100-120 V	900033893
102A	Circuit card Mains input 220-240 V	900033885
102B	Circuit card Mains input 100-120 V	900033892
103	Support	900033888

Position	Product name	Item number
104	Temperature sensor	900033883
105A	Heating plate 220-240 V	900033881
105B	Heating plate 100-120 V	900033891
106A	Ultrasound well 220-240 V	900033880
106B	Ultrasound well 100-120 V	900033890
107	Fuse 5x20mm, T15A	900033837
108	Equipment base	900033836
501	Plastic lid	18500001
502	Stainless steel insert cover	18500002
503	Stainless steel basket	18500003
504	Stainless steel immersion basket 59mm	18500004
505	Plastic bath insert for acids	18500005
506	Cleaning jar	18500006
507	Plastic container with lid	18500007 *

\* Wear and tear parts, excluded from the guarantee.