


	Upper anteriors	Lower anteriors	Posteriors	0°	20°
○	● O 20	L 30		30, 32	32
	● O 30	L 40		32	32
◻	■ R 20	L 20		32	32, 34
	■ R 30	L 25		32	36
	■ R 40	L 55		32	36
◻	■ S 05	L 10		30	32
	■ S 10	L 10		30	32
	■ S 20	L 20		30, 32	32
	■ S 30	L 40		32	32, 34
	■ S 35	L 55		32	36
	■ S 40	L 55		32	36
△	▲ T 10	L 10		30	32
	▲ T 20	L 30		30, 32	32
	▲ T 30	L 20		32	32
	▲ T 35	L 40		32	32
	▲ T 40	L 50		32	34
◻	▲ TO 05	L 10		30	32
	▲ TO 10	L 15		30	32
	▲ TO 20	L 20		30, 32	32
	▲ TO 30	L 40		32	32
	▲ TO 40	L 45		32	34



This mould chart complies with EN ISO 22112:2017. / Diese Formenkarte erfüllt EN ISO 22112:2017. / Cette carte forme remplit EN ISO 22112:2017. / Esta carta de formas es conforme con la EN ISO 22112:2017. / Questa tabella forme soddisfa EN ISO 22112:2017. / Esta carta molde está em conformidade com EN ISO 22112:2017. / Deze vormenkaart voldoet aan EN ISO 22112:2017. / Detta formkort uppfyller EN ISO 22112:2017. / Dette formkort opfylder EN ISO 22112:2017. / Tämä muotokortti täyttää EN ISO 22112:2017. / Dette formkortet oppfyller EN ISO 22112:2017. / See vormikaart vastab standardile EN ISO 22112:2017. / Ši formas karte atbilst EN ISO 22112:2017. / Ši formas korta išpildo EN ISO 22112:2017. / Ta karta ksztaftu spełnia EN ISO 22112:2017. / Þetta mótakort er í samræmi við EN ISO 22112:2017. / Ova kartica oblika ispunjava EN ISO 22112:2017. / Тази матрица е в съответствие с EN ISO 22112:2017. / Tato karta tvaru splňuje EN ISO 22112:2017. / Tato tvarová karta splňa EN ISO 22112:2017. / Ова картица калупа је у складу са EN ISO 22112:2017. / Оваа карта за мувла е во согласност со EN ISO 22112:2017. / 이 몰드 카드는 EN ISO 22112:2017을 준수합니다. / يتوافق مخطط القوالب هذا مع EN ISO 22112:2017.



Manufacturer

Kulzer GmbH
Leipziger Straße 2
63450 Hanau
Germany

Contact in Germany

Kulzer GmbH
Leipziger Straße 2
63450 Hanau
Germany
info.lab@kulzer-dental.com

artic digital

- GB Mould chart
- DE Formenkarte
- FR Carte des Formes
- ES Carta de Formas
- IT Tabella forme
- PT Carta Molde
- NL Vormenkaart
- SE Formkort
- DK Formkort
- FI Muotto korti
- NO Formkort
- EE Vormikaart
- LV Formas karte
- LT Formos kortelė
- PL Karta kształtu
- IS Myglatöflu
- HR Kartica oblika
- BG плесен диаграма
- CZ Karta tvaru
- SK Tvarová karta
- RS картица калупа
- MK форма картичка
- KR 모양 카드
- ARA مخطط القالب



Artic® Digital





































Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP













artic digital 6

Anteriors / Frontzähne / Antérieures / Anteriores / Anteriori / Anteriores / Tandem /
 Framtänder / For tænder / Etohampaat / Fortennene / Eesmised hambad /
 Priekšējie zobi / Priekiniai dantys / Zęby przednie / Framhliðar / Prednji zubi / предни
 зъби / Přední zuby / Predné zuby / предњи зуби / предните заби / 앞니 / الأمامية

O 20			46,1 9,7 8,4	S 35			50,6 10,8 9,2
O 30			46,6 10,9 8,4	S 40			52,3 12,4 9,8
R 20			42,9 9,9 7,8	T 10			44,1 9,0 8,2
R 30			47,7 12,1 8,8	T 20			45,6 9,8 8,8
R 40			51,1 13,3 9,5	T 30			45,0 10,8 8,1
S 05			40,3 8,9 7,4	T 35			47,4 10,0 8,8
S 10			44,2 9,1 7,9	T 40			49,1 12,1 8,7
S 20			44,0 9,9 7,9	TO 05			42,1 9,0 7,8
S 30			48,2 10,9 8,6	TO 10			43,3 9,3 8,3





















artic digital 6

Anteriors / Frontzähne / Antérieures / Anteriores / Anteriori / Anteriores / Tandem /
 Framtänder / For tænder / Etohampaat / Fortennene / Eesmised hambad /
 Priekšējie zobi / Priekiniai dantys / Zęby przednie / Framhliðar / Prednji zubi / предни
 зъби / Přední zuby / Predné zuby / предњи зуби / предните заби / 앞니 / الأمامية

TO 20		43,6 10,2 7,7	L 25		36,8 9,3 5,5
TO 30		45,9 10,4 8,6	L 30		35,9 9,0 5,5
TO 40		47,8 11,0 8,9	L 40		36,0 9,0 5,5
L 10		33,5 7,8 5,1	L 45		38,9 9,6 5,9
L 15		33,9 9,5 5,2	L 50		39,2 11,1 5,9
L 20		35,2 9,1 5,4	L 55		40,2 10,6 6,0

artic digital 8 0°/20°

Posteriors/Seitenzähne/Postérieures/Posteriores/Posteriori/Posteriores/Kiezen/
 Bakre tänder/Bageste tænder/Takimmaisest hampaat/Bakre tenner/Tagumised
 hambad/Aizmugurējie zobi/Užpakaliniai dantys/Zęby tylne/Bakhliðar/Stražnji zubi/
 задни зъби/Zadní zuby/Zadné zuby/задњи зуби/задните заби/구치부/الخلفية

30 U	0°			27,2 7,9
30 L				28,9 8,6
32 U	0°			33,1 8,8
32 L				34,8 9,2
32 U	20°			31,1 9,1
32 L				33,9 8,0
34 U	20°			32,9 9,4
34 L				35,7 8,3
36 U	20°			34,7 10,0
36 L				37,4 9,0